



trio  smart[®]
SAMPLE COLLECTION INSTRUCTIONS

WATCH THE SAMPLE COLLECTION
INSTRUCTIONS VIDEO

triosmartbreath.com/instructions



TABLE OF CONTENTS

WHAT'S IN YOUR KIT	3
COLLECT YOUR SUGAR SUBSTRATE DRINK	5
PREPARE FOR YOUR SAMPLE COLLECTION	6
COMPLETE YOUR REQUISITION FORM	7
COLLECT YOUR BREATH SAMPLES	8
SHIP YOUR SAMPLES FOR TESTING	15
CONTACT INFORMATION	16

WELCOME TO **trio** **smart**®

The **trio-smart** breath test measures the levels of hydrogen, methane, and hydrogen sulfide in a patient's breath after the consumption of a sugar substrate drink.

THE COMPLETE FERMENTED GAS PROFILE OF THE MICROBIOME

HYDROGEN

Indicative of:
**Small Intestinal Bacterial
Overgrowth (SIBO)**

Correlated with:
No correlation with symptoms

METHANE

Indicative of:
**Intestinal Methanogenic
Overgrowth (IMO)**

Correlated with:
Constipation



HYDROGEN SULFIDE

Indicative of:
Excess Hydrogen Sulfide

Correlated with:
Diarrhea


WHAT'S IN YOUR KIT?



One biohazard bag
for the safe return
of your breath
samples



One mouthpiece
that attaches to
the valves on the
bags



**Nine breath
collection bags**
with white valves



One requisition form
to be completed
and returned
with your kit



Nine sticker labels
for you to fill out
and attach
to the bags



**One pre-paid
shipping label**
for returning
the kit

COLLECT YOUR SUGAR SUBSTRATE DRINK

Your doctor has selected the type of substrate drink for you to consume during sample collection.



You may have received **a signed prescription form for lactulose** in your trio-smart kit, or your doctor may have sent the prescription directly to your pharmacist. Please fill the prescription prior to taking the test. You are required to drink a single dose of 10gm of lactulose in a 15ml solution (approximately 1 tablespoon) with or followed by one cup of water **after collecting your first sample**. Discard any excess solution.

If you received a bottle of glucose in your kit, please continue to the next step. You are required to drink the full dose of 75gm of glucose solution **after collecting your first sample**.

PREPARE FOR YOUR SAMPLE COLLECTION

The North American Consensus on Breath Testing for GI Disorders recommends the following preparation before sample collection.



If tolerated by the patient, **promotility drugs and laxatives should be stopped** at least 7 days prior to sample collection.



A restricted diet should be followed on the day prior to sample collection.

▶ During this time period, restrict your diet to only meat (poultry, fish, beef or pork), tofu, eggs, potatoes and white rice, limiting your seasonings to oil, salt and pepper. In addition to water, drink only coffee or tea without any additives. This diet has been specially designed to ensure the accuracy of your test results. **Do not eat or drink anything other than the items listed here.**



A **fasting period of 8–12 hours** should be observed prior to sample collection.

▶ Do not eat during sample collection.

COMPLETE YOUR REQUISITION FORM

- 1 Make sure the **Patient Information** and **Patient's Insurance Information**, including the requested **policy holder information** portions of the enclosed requisition form, are complete and accurate.
- 2 Please attach **copies of the front and back of your insurance card**. If you prefer, you can email images of your card to support@triosmartbreath.com.
- 3 Sign and date the requisition form on the line marked **Patient Signature**.




Remember to enclose the requisition form and copies of your insurance card with your samples when shipping them to the lab.

COLLECT YOUR BREATH SAMPLES

STEP 1: PREPARE AND APPLY THE STICKER LABEL.




Complete a sticker label with your **sample number** (1, 2, 3, etc.), **sample collection time**, **date of collection**, **name**, and **date of birth**.



Attach the label to a deflated sample collection bag.


STEP 2: ATTACH THE MOUTHPIECE AND OPEN THE VALVE.



Attach the mouthpiece to the end of the valve until it fits snugly. The **small hole** on the mouthpiece should be **facing up**.

Turn the cap on the valve **counterclockwise** until it is fully in the “**Open**” position.


STEP 3: COLLECT YOUR BREATH SAMPLE.



Take a deep breath and **exhale continuously** into the mouthpiece. Air will escape from the hole on the top of the mouthpiece. **After two seconds**, while continuing to exhale, **cover the hole** with your finger. The bag will then inflate.

► Only exhale **once** into each collection bag.


STEP 4: CLOSE THE VALVE AND REMOVE THE MOUTHPIECE.



Close the cap by turning it **clockwise** until it is fully in the **“Closed”** position.

Remove the mouthpiece from the valve.

STEP 5: CONSUME THE SUGAR SUBSTRATE DRINK AND WAIT 15 MINUTES.




After collecting your first breath sample, please drink your prescribed sugar substrate drink (lactulose or glucose). **Wait 15 minutes** before collecting your next sample.

► If prescribed lactulose, please drink a single dose of 10gm of lactulose in a 15ml solution (approximately 1 tablespoon) with or followed by one cup of water. Discard any excess solution.

► If prescribed glucose, please drink the full 75gm dose of glucose solution provided in your kit.

STEP 6: COLLECT YOUR REMAINING SAMPLES.

Repeat steps 1–4 for each remaining bag to collect eight more breath samples. Remember to **wait 15 minutes** between the collection of each breath sample.



► Confirm that each label is completed with the **time of collection** and is properly applied to each bag.

SHIP YOUR BREATH SAMPLES FOR TESTING



Place all filled sample collection bags into the provided biohazard bag and seal it tightly.



Place the biohazard bag, along with your signed requisition form and copies of your insurance card, inside of the trio-smart box.



Place the pre-paid return shipping label over the original shipping label, and close the box with packing tape.



Drop the kit off at any UPS drop-off location ([ups.com/dropoff](https://www.ups.com/dropoff)) on any day of the week, no later than 24 hours after sample collection.

▶ Results will be available within **seven business days** after your sample is received.

QUESTIONS? COMMENTS? FEEDBACK?

Website:

triosmartbreath.com

Email:

support@triosmartbreath.com

Phone:

(855) 888-1230

trio  smart[®]

